

Product Information

MemDX™ Membrane Protein Human MCEMP1 (Mast cell expressed membrane protein 1)

Full Length

Cat. No.: **MPC2781K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MCEMP1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MCEMP1

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEVEEIIYKHQEVKMQAPAFRDKKQGVSAKNQGAHDPDYENITLAFKNQDH
AKGGHSRPTSQVPAQCRPPSDSTQVPCWLYRAILSLYILLALAFVLCIIL
SAFIMVKNAEMSKELLGFKRELWNVSNSVQACEERQKRGWDSVQQSITMV
RSKIDRLLETTLAGIKNIDTKVQKILEVLQKMPQSSPQ

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MCEMP1

Full Name

Mast cell expressed membrane protein 1

Introduction

This gene encodes a single-pass transmembrane protein. Based on its expression pattern, it is speculated to be involved in regulating mast cell differentiation or immune responses.

Alternative Names

MCEMP1; C19orf59; epididymis secretory sperm binding protein; Mast cell expressed membrane protein 1

Gene ID

[199675](#)

UniProt ID

[Q8IX19](#)