

Product Information

MemDX™ Membrane Protein Human STOML3 (Stomatin like 3) Full Length

Cat. No.: **MPC2782K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human STOML3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

STOML3

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MDSRVSSPEKQDKENFVGVNNKRLGVCGWILFSLSFLLVIITFPISIWMC
LKIIEYERAVVFR LGRIQADKAKGPGLILVLP CIDVFVKVDLRTVTCNI
PPQEILTRDSVTTQVDGVVYYRIYSAVSAVANVNDVHQATFLLAQTTLRN
VLGTQTLSQLAGREEIAHSIQTLDDATELWGIRVARVEIKDVRIPVQL
QRSMAAAEAEATREARAKVLAAEGEMNASKSLKSASMVLAESPIALQLRYL
QTLSTVATEKNSTIVFPLPMNILEGIGGVSYDNHKKLPNKA

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

STOML3

Full Name

Stomatin like 3

Alternative Names

STOML3; SRO; Epb7.2I; stomatin-like protein 3; EPB72-like 3; SLP-3; erythrocyte band 7 integral membrane protein; protein 7.2B; stomatin (EPB72)-like 3; Stomatin like 3

Gene ID

[161003](#)

UniProt ID

[Q8TAV4](#)