

Product Information

MemDX™ Membrane Protein Human AJAP1 (Adherens junctions associated protein 1) Full Length

Cat. No.: **MPC2802K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human AJAP1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

AJAP1

Protein Length

Full length

Protein Class

Cell adhesion

TMD

1

Sequence

MWIIQQLLGLSSMSIRWPGRPLGSHAWILIAMFQLAVDLPACEALGPGPEF
WLLPRSPRPPRLWSFRSGQPARVPAPVWSPRPPRVERIHGQMMPRARR
AHRPRDQAAALVPKAGLAKPPAAAKSSPSLASSSSSSSSAVAGGAPEQQA
LLRRGKRHLQGDGLSSFDSRGSRPTTETEFIAWGPTGDEEALESNTFPGV
YGPTTVSILQTRKTTVAATTTTTTTATPMTLQTKGFESLDPRRRRIPGGV
STTEPSTSPSNNGEVTQPPRILGEASGLAVHQIITITVSLIMVIAALITT
LVLKNCCAQSGNTRRNHQKTNQQEESCQNLTDFPSARVPSSLDIFTAY
NETLQCSHECVRASVPVYTDETLHSTTGEYKSTFNNGRNPSSSDRHLIPVA
FVSEKWFEISC

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

AJAP1

Full Name

Adherens junctions associated protein 1

Alternative Names

AJAP1; MOT8; SHREW1; SHREW-1; adherens junction-associated protein 1; membrane protein shrew-1; transmembrane protein SHREW1; Adherens junctions associated protein 1

Gene ID

[55966](#)

UniProt ID

[Q9UKB5](#)