

Product Information

MemDX™ Membrane Protein Human SLC35F1 (Solute carrier family 35 member F1) Full

Length

Cat. No.: **MPC2812K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SLC35F1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SLC35F1

Protein Length

Full length

Protein Class

Transporter

TMD

10

Sequence

MIPPEQPQQQLQPPSPAPPNHVVTIENLPAEGSGGGGSLSSASSRAGVRQ
RIRKVLNREMLISVALGQVLSLLICGIGLTSKYLSEDFHANTPVFQSFLN
YILLFLVYTTTTLAVRQGEENLLAILRRRWWKYMILGLIDLEANLVVKAY
QYTTLTSLQLLDCFVIPVVILLSWFFLLIRYKAVHFIGIVVCILGMGCMV
GADVLVGRHQGAGENKLVGDLLVLGGATLYGISNVWEEYIIRTLRVEFL
GMIGLFGAFFSGIQLAIMEHKELLKVPWDWQIGLLYVGFSACMFGLYSFM
PVVIKKTSAATSVNLSLLTADLYSLFCGLFLFHYKFSGLYLLSFFTILIGL
VLYSSTSTYIAQDPRVYKQFRNPSGPVVDLPTTAQVEPSVTYTSLGQETE
EEPHVRVA

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SLC35F1

Full Name

Solute carrier family 35 member F1

Alternative Names

SLC35F1; C6orf169; dJ230I3.1; Solute carrier family 35 member F1

Gene ID

[222553](#)

UniProt ID

[Q5T1Q4](#)