

Product Information

MemDX™ Membrane Protein Human CD72 (CD72 molecule) Full Length

Cat. No.: **MPC2845K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CD72 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD72

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MAEAITYADLRFVKAPLKSSISRLGQDPGADDDGEITYENVQVPAVLGV
PSSLASSVLGDKAIVKSEQPTASWRAVTSPAVGRILPCRTTCLRYLLLGL
LLTCLLLGVTAICLGVRYLQVSQQQLQQTNRVLEVTNSSLRQQRLRLKITQL
GQSAEDLQGSRRRLAQSQEALQVEQRAHQAAEGQLQACQADRQKTKETLQ
SEEQRRRALEQKLSNMENRLKPFFTCGSADTCCPSGWIMHQKSCFYISLT
SKNWQESQKQCETLSSKLATFSEIYPQSHSYFLNSLLPNGGSGNSYWTG
LSSNKDWKLTDDTQRTRTYAQSSKCNKVHKTWSWWTLESESCRSSLPYIC
EMTAFRFPD

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

CD72

Full Name

CD72 molecule

Alternative Names

CD72; LYB2; CD72b; B-cell differentiation antigen CD72; CD72 antigen; lyb-2; CD72 molecule

Gene ID

[971](#)

UniProt ID

[P21854](#)