

Product Information

MemDX™ Membrane Protein Human PRAF2 (PRA1 domain family member 2) Full Length

Cat. No.: **MPC2856K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human PRAF2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

PRAF2

Protein Length

Full length

Protein Class

Transporter

TMD

4

Sequence

MSEVRLPPLRALDDFVLGSARLAAPDCDPQRWCHRVINNLLYYQTNYLL
CFGIGLALAGYVRPLHTLLSALVVAVALGVLVWAAETRAAVRRCRRSHPA
ACLA AVLAVGLLVLVWAGGACTFLFSIAGPVLLILVHASLRLRLNKNKIE
NKIESIGLKRTMPMGLLLEALGQEQEAGS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

PRAF2

Full Name

PRA1 domain family member 2

Alternative Names

PRAF2; JM4; Yip6a; PRA1 family protein 2; Jena-Muenchen 4; PRA1 domain family member 2

Gene ID

[11230](#)

UniProt ID

[O60831](#)