

Product Information

MemDX™ Membrane Protein Human HBEGF (Heparin binding EGF like growth factor) Full Length

Cat. No.: **MPC2937K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human HBEGF membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

HBEGF

Protein Length

Full length

Protein Class

Growth factor

TMD

1

Sequence

MKLLPSVVLKFLAAVLSALVTGESLERLRRGLAAGTSNPDPTVSTDQL
LPLGGGRDRKVRDLQEADLDLLRVTLSKPKALATPNKEEHGKRKKKGKG
LGKKRDPCLRKYKDFCIHGECKYVKELRAPSCICHPGYHGERCHGLSLPV
ENRLTYDHTTILAVVAVVLSSVCLLVIVGLLMFRYHRRGGYDVENEKV
KLGMTNSH

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

HBEGF

Full Name

Heparin binding EGF like growth factor

Alternative Names

HBEGF; DTR; DTS; DTSF; HEGFL; proheparin-binding EGF-like growth factor; diphtheria toxin receptor (heparin-binding EGF-like growth factor); diphtheria toxin receptor (heparin-binding epidermal growth factor-like growth factor); heparin-binding epidermal growth factor; Heparin binding EGF like growth factor

Gene ID

[1839](#)

UniProt ID

[Q99075](#)