

Product Information

MemDX™ Membrane Protein Human CD1A (CD1a molecule) Full Length

Cat. No.: **MPC2941K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CD1A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD1A

Protein Length

Full length

Protein Class

Immunity

TMD

1

Sequence

MLFLLLPLLA VLP GDGNADGLKEPLSFHVTWIASFYNH SWKQNLVSGWLS
DLQTHTWDSNSSTIVFLCPWSRGNF SNEEWKELET LFRIRTIRSFEGIRR
YAH ELQFEY PFEIQVTGGCE LHSGKVSGSFLQLAYQGSDFVSFQNN SWLP
YPVAGNMAKH FCKVLNQNHENDITHNLLSDTCPRFILGLLDAGKAHLQR
QVKPEAWLSHG PPGPGHLQLVCHVSGFY PKPVWVMW MRGEQEQQGTQRG
DILPSADGTWYLRATLEVAAGEAADLSCR VKHSSLEGQDIVLYWEHHSSV
GFILAVIVPLLLLIGLALWFRKRCFC

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

CD1A

Full Name

CD1a molecule

Introduction

This gene encodes a member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene localizes to the plasma membrane and to recycling vesicles of the early endocytic system. Alternative splicing results in multiple transcript variants.

Alternative Names

CD1A; R4; T6; CD1; FCB6; HTA1; T-cell surface glycoprotein CD1a; CD1A antigen, a polypeptide; T-cell surface antigen T6/Leu-6; cluster of differentiation 1 A; cortical thymocyte antigen CD1A; differentiation antigen CD1-alpha-3; epidermal dendritic cell marker CD1a; hTa1 thymocyte antigen; CD1a molecule

Gene ID

[909](#)

UniProt ID

[P06126](#)