

Product Information

MemDX™ Membrane Protein Human SPCS3 (Signal peptidase complex subunit 3) Full

Length

Cat. No.: **MPC3002K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SPCS3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SPCS3

Protein Length

Full length

Protein Class

Protease

TMD

1

Sequence

MNTVLSRANSLFAFSLSVMAALTFGCFITAFKDRSVPVRLHVSRIMLKN
VEDFTGPRERSDLGFITFDITADLENIFDWNVKQLFLYLSAEYSTKNNAL
NQVVLWDKIVLRGDNPKLLLKDMKTKYFFFDDGNGLGKGNRNVTLTSLWNV
VPNAGILPLVTGSGHVSVPFPDYEITKSY

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SPCS3

Full Name

Signal peptidase complex subunit 3

Alternative Names

SPCS3; SPC3; SPC22; SPC23; PRO3567; YLR066W; SPC22/23; SPase 22 kDa subunit; SPase 22/23 kDa subunit; microsomal signal peptidase 22/23 kDa subunit; microsomal signal peptidase 23 kDa subunit; signal peptidase complex subunit 3 homolog; Signal peptidase complex subunit 3

Gene ID

[60559](#)

UniProt ID

[P61009](#)