

Product Information

MemDX™ Membrane Protein Human CD1C (CD1c molecule) Full Length

Cat. No.: **MPC3016K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CD1C membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD1C

Protein Length

Full length

Protein Class

Immunity

TMD

1

Sequence

MLFLQLLLALLLPGGDNADASQEHVSFHVLIQIFS FVNQSWARGQGSGWL
DELQTHGWDSSEGTIIFLHNWSKGNFSNEELSDLELLFRFYLFGLTREIQ
DHASQDYSKYPFEVQVKAGCELSGKSPEGFFQVAFNGLDLLSFQNTTWV
PSPGCGSLAQSVCHLLNHQYEGVTETVYNLIRSTCPRFLLGLLDAGKMYV
HRQVRPEAWLSSRPSLGSGQLLVCHASGFYKPKVWVTWMRNEQEQLGTK
HGDILPNADGTWYLQVILEVASEEPAGLSCRVRHSSLGGQDIILYWGHHF
SMNWIALVVIVPLVILIVLVLWFKKHCSYQDIL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

CD1C

Full Name

CD1c molecule

Introduction

This gene encodes a member of the CD1 family of transmembrane glycoproteins, which are structurally related to the major histocompatibility complex (MHC) proteins and form heterodimers with beta-2-microglobulin. The CD1 proteins mediate the presentation of primarily lipid and glycolipid antigens of self or microbial origin to T cells. The human genome contains five CD1 family genes organized in a cluster on chromosome 1. The CD1 family members are thought to differ in their cellular localization and specificity for particular lipid ligands. The protein encoded by this gene is broadly distributed throughout the endocytic system via a tyrosine-based motif in the cytoplasmic tail. Alternatively spliced transcript variants of this gene have been observed, but their full-length nature is not known.

Alternative Names

CD1C; R7; CD1; CD1A; BDCA1; T-cell surface glycoprotein CD1c; CD1C antigen, c polypeptide; cortical thymocyte antigen CD1C; differentiation antigen CD1-alpha-3; CD1c molecule

Gene ID

[911](#)

UniProt ID

[P29017](#)