

Product Information

MemDX™ Membrane Protein Human ZPLD1 (Zona pellucida like domain containing 1) Full Length

Cat. No.: **MPC3444K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZPLD1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZPLD1

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEQIWLLLLLTIRVLPGSAQFNGYNCDANLHSRFP AERDISVYCGVQAIT
MKINFCTVLFSGYSETDLALNGRHGDSHCRGFINNNTFPAVVIFIINLST
LEGCGNNLVVSTIPGVSAYGNATSVQVGNISGYIDTPDPPTIISYLPGLL
YKFSCSYPLEYLVNNTQLASSSAAISVRENNGTFVSTLNLLLYNDSTYNQ
QLIIPSIGLPLKTKVFAAVQATNLDGRWNVLMDYCYTTPSGNPND DIRYD
LFLSCDKDPQTTVIENGRSQRGRFSFEVFRFVKHKNQKMSTVFLHCVTKL
CRADDCPFLMPICSHRERRDAGRRTTWSPQSSSGSAVLSAGPIITRSD ET
PTNNSQLGSPSPMPFQLNAITSALISGMVILGVTSFSLLLCSLALLHRKG
PTSLVLNGIRNPVFD

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZPLD1

Full Name

Zona pellucida like domain containing 1

Alternative Names

ZPLD1; zona pellucida-like domain-containing protein 1; ZP domain-containing protein 1; cupulin; Zona pellucida like domain containing 1

Gene ID

[131368](#)

UniProt ID

[Q8TCW7](#)