

Product Information

MemDX™ Membrane Protein Human TRGV9 (T cell receptor gamma variable 9) Full Length

Cat. No.: **MPC3446K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TRGV9 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TRGV9

Protein Length

Full length

Protein Class

Receptor; Immunity

Sequence

MLSLLHTSTLAVLALCVYAGHLEQPQISSTKLSKTARLECVSGITI
SATSVYWYRERPGEVIQFLVSIYDGTVRKESGIPSGKFEVD RIPETSTS
TLTIHNVEKQDIATYYCALWEV

Product Description

Expression Systems

E.coli expression system

Tag

Based on specific requirements

Protein Format

Soluble

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TRGV9

Full Name

T cell receptor gamma variable 9

Alternative Names

TRGV9; V2; TRGC1; TCRGV9; T-cell receptor, gamma, variable region V9; T cell receptor gamma variable 9

Gene ID

[6983](#)

UniProt ID

[Q99603](#)