

Product Information

MemDX™ Membrane Protein Human TRGC1 (T cell receptor gamma constant 1) Full Length

Cat. No.: **MPC3450K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TRGC1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TRGC1

Protein Length

Full length

Protein Class

Receptor; Immunity

TMD

1

Sequence

DKQLDADVSPKPTIFLPSIAETKLQKAGTYLCLLEKFFPDVIKHWQEKK
SNTILGSQEGNTMKTNDTYMKFSWLTVPKSLDKEHRCIVRHENNKNGVD
QEIIFFPIKTDVITMDPKDNCSDANDTLLLQLTNTSAYMYLLLLLKSV
VYFAITCCLLRRTAFCCNGEKS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

TRGC1

Full Name

T cell receptor gamma constant 1

Alternative Names

TRGC1; C1; TCRG; TCRGC1; T-cell receptor, gamma, constant region C1; T cell receptor gamma constant 1

Gene ID

[6966](#)

UniProt ID

[P0CF51](#)