

Product Information

MemDX™ Membrane Protein Human GPR32P1 (GPR32 pseudogene 1) Full Length

Cat. No.: **MPC3454K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human GPR32P1 membrane protein expressed in Baculovirus/Insect expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR32P1

Protein Length

Full length

Protein Class

GPCR

TMD

5

Sequence

MNGVSEGTRGCSDRQPGALTQGHSCSRKMNASRCLSEEVGSLRPLTMAVL
SASFVVGVLGNGLVPWVTVFRMARTVSTVCFFHLALADFMLSLSLPILVY
YIVSRQWLLGEWACKLYTGfVFLTfSTSNCLLVISVDRCISVLYPVWAL
NHRTEQRASWLAFGVWLLAAALCSAHLKFRTRKWNGCMQCYLQFNLENE
TAQMWTQEVFGRQMAVIMAHFLLGFLGPLAIIGTCAHLIRAKLLREGWVH
ANRPKRLLLLVLVSALSAGSHLT

Product Description

Expression Systems

Baculovirus/Insect expression system

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR32P1

Full Name

GPR32 pseudogene 1

Alternative Names

GPR32P1; GPR32P; hGPCR39; G protein-coupled receptor 32, pseudogene 1; Putative G-protein coupled receptor GPCR39; GPR32 pseudogene 1

Gene ID

[2855](#)

UniProt ID

[Q8NGA4](#)