

# Product Information

## MemDX™ Membrane Protein Human MS4A13 (Membrane spanning 4-domains A13) Full

### Length

Cat. No.: **MPC3476K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MS4A13 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

MS4A13

#### Protein Length

Full length

#### Protein Class

Transporter

#### TMD

4

#### Sequence

MIGIFHIFMWYFLLVLYMGQIKGAFGTYEPVYKTGCTLWGIFFIAGVF  
LIRVTKYPTRSGIISTLIINIICITITAVTLTIELSHFNSVSYRNYG  
QAKLGREVSRILLFFYGLEFSIALTHSIYSCSNLFRQNDLTSVTEEAES  
TP

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

#### Form

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

MS4A13

**Full Name**

Membrane spanning 4-domains A13

**Alternative Names**

MS4A13; membrane-spanning 4-domains subfamily A member 13; testis-expressed transmembrane protein 4; testis-expressed transmembrane-4 protein; Membrane spanning 4-domains A13

**Gene ID**

[503497](#)

**UniProt ID**

[Q5J8X5](#)