

Product Information

MemDX™ Membrane Protein Human ZDHHC23 (Zinc finger DHHC-type palmitoyltransferase

23) Full Length

Cat. No.: **MPC3480K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC23 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC23

Protein Length

Full length

Protein Class

Transferase

TMD

7

Sequence

MTQKGSMKPVKKKKTEEPELEPLCCCEYIDRNGEKNHVATCLCDCQDLDE
GCDRWITCKSLQPETCERIMDTISDRLRIPWLRGAKKVNISIIPPLVLLP
VFLHVASWHFLLGVVLTSLPVLALWYYLTHRRKEQTLFFLSLGLFSLG
YMYVFLQEVVPKGRVGPVQLAVLTCGLFLILLALHRAKKNPGYLSNPAS
GDRSLSSSQLECLSRKGQEKTKGFPADMSGSLNNRTTKDDPKGSSKMPA
GSPTKAKEDWCAKCQLVRPARAWHCRICGICVRRMDHHCWVWINSVGESN
HQAIFILALLIFLLTSVYGITLTLDITCRDRSVFTALFYCPGVYANYSSAL
SFTCVWYSVIITAGMAYIFLIQLINISYNVTEREVQQALRQKTGRRLLCG
LIVDTGLLG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC23

Full Name

Zinc finger DHHC-type palmitoyltransferase 23

Alternative Names

ZDHHC23; NIDD; DHHC-23; palmitoyltransferase ZDHHC23; zinc finger DHHC domain-containing protein 23; zinc finger DHHC-type containing 23; Zinc finger DHHC-type palmitoyltransferase 23

Gene ID

[254887](#)

UniProt ID

[Q8IYP9](#)