

## Product Information

### MemDX™ Membrane Protein Human MEGF9 (Multiple EGF like domains 9) Full Length

Cat. No.: **MPC3550K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MEGF9 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

#### Product Specifications

##### Host Species

Human

##### Target Protein

MEGF9

##### Protein Length

Full length

##### Protein Class

Receptor

##### TMD

1

##### Sequence

MNGGAERAMRSLPSLGGLALLCCAAAAAASAGNVTGGGGAAGQ  
VDASPGPGLRGEP SHPFRATAPTAQAPRTGPPRATVHRPLAATSPAQSP  
ETTPLWATAGPSSSTTFQAPLGPSPTTPPAAERTSTTSQAPTRPAPTTLST  
TTGPAPTPPVATTVPAPTTPTPTDLPSSSNSSVLPTPPATEAPSSPPP  
EYVCNCSVVGSLNVNRCNQTTGQCECRPGYQGLHCETCKEGFYLNYSGL  
CQPCDCSPHGALSIPCNSSGKCQCKVGVIGSICDRCQDGYGFSKNGCLP  
CQCNNRSASCDALTGACLNQCENSKGNHCEECKEGFYQSPDATKECLRCP  
CSAVTSTGSCSIKSELEPECDQCKDGYIGPNCNKCENGYYNFDSICRKC  
QCHGHVDPVKT PKICKPESGECINCLHNTTGFWCENCLEGYVHDLEGNCI  
KKEVILPTPEGSTILVSNASLTTSVPTPVINSTFTPTTLQTIFSVSTSEN  
STSALADVSWTQFNIIILTVIIIVVLLMGFVGAVMYREYQNRKLNAPF  
WTIELKEDNISFSSYHDSIPNADVSGLLEDDGNEVAPNGQLTLTPIHNY  
KA

#### Product Description

##### Expression Systems

HEK293

##### Tag

Based on specific requirements

**Protein Format**

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

MEGF9

**Full Name**

Multiple EGF like domains 9

**Alternative Names**

MEGF9; EGFL5; multiple epidermal growth factor-like domains protein 9; EGF-like-domain, multiple 5; epidermal growth factor-like protein 5; Multiple EGF like domains 9

**Gene ID**

[1955](#)

**UniProt ID**

[Q9H1U4](#)