

Product Information

MemDX™ Membrane Protein Human AMIGO3 (Adhesion molecule with Ig like domain 3) Full Length

Cat. No.: **MPC3551K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human AMIGO3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

AMIGO3

Protein Length

Full length

Protein Class

Cell adhesion

TMD

1

Sequence

MTWLVLLGTLLCMLRVGLGTPDSEGFPPRALHNCYPYKCICAADLLSCTGL
GLQDVPALPAATADLDLSHNAQLRPGWLAPLFQLRALHLDHNELDAL
GRGVFVNASGLRLLDLSSNLRALGRHLDGLGALEKLLLFNNRLVHLDE
HAFHGLRALSHLYLGCNELASFSDHLHGLSATHLLTDLSSNRLGHISV
PELAALPAFLKNGLYLHNNPLPCDRLYHLLQRWHQRGLSAVRDFAREYV
CLAFKVPASRVRFQHSRVFENCSSAPALGLERPEEHLYALVGRSLRLYC
NTSVPAMRIAWVSPQQELLRAPGSRDGSIAVLADGSLAIGNVQEQHAGLF
VCLATGPRLLHNNQTHEYNVSVHFPRPEPEAFNTGFTTLLGCAVGLVLVLL
YLFAPPCRCRCRRACRCRRWPQTSPQLQELSAQSSVLSTTPPDAPSRKASV
HKHVVFLEPGRRGLNGRVQLAVAEFDLYNPGGLQLKAGSESASSIGSEG
PMTT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

AMIGO3

Full Name

Adhesion molecule with Ig like domain 3

Alternative Names

AMIGO3; ALI3; AMIGO-3; amphoterin-induced protein 3; alivin-3; amphoterin-induced gene and ORF 3; Adhesion molecule with Ig like domain 3

Gene ID

[386724](#)

UniProt ID

[Q86WK7](#)