

Product Information

MemDX™ Membrane Protein Human FXYP6P3 (FXYP domain containing ion transport regulator 6 pseudogene 3) Full Length

Cat. No.: **MPC3605K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human FXYP6P3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FXYP6P3

Protein Length

Full length

Protein Class

Transporter; Ion channel

TMD

1

Sequence

MEVVLIFVYSLVPVVLASAAKEKEIDPFHYNQTLRIGGLVFDVVLFLV
PSCHLLSHRCKCSFNQKPQDPGDKEAQVENFITANAKEPQKAKN

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

FXYP6P3

Full Name

FXYP domain containing ion transport regulator 6 pseudogene 3

Alternative Names

FXYP6P3; FXYP8; FXYP domain containing ion transport regulator 8; FXYP domain containing ion transport regulator 6 pseudogene 3

Gene ID

[406875](#)

UniProt ID

[P58550](#)