

# Product Information

## MemDX™ Membrane Protein Human RABAC1 (Rab acceptor 1) Full Length

Cat. No.: **MPC3615K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human RABAC1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

RABAC1

#### Protein Length

Full length

#### Protein Class

Transporter

#### TMD

4

#### Sequence

MAAQKDQQKDAEAEGLSGTTLLPKLIPSGAGREWLERRRATIRPWSTFVD  
QQRFSRPRNLGELCQRLVRNVEYYQSNYVFVFLGLILYCVVTSPMLLVAL  
AVFFGACYILYLRTLESKLVLFGREVSPAHHQYALAGGISFPFFWLAGAGS  
AVFWVLGATLVVIGSHAAFHQIEAVDGEELQMEPV

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

#### Form

Liquid

#### Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

## Target

### Target Protein

RABAC1

### Full Name

Rab acceptor 1

### Alternative Names

RABAC1; PRA1; YIP3; PRAF1; prenylated Rab acceptor protein 1; PRA1 domain family 1; PRA1 family protein 1; Rab acceptor 1 (prenylated); prenylated Rab acceptor 1; Rab acceptor 1

### Gene ID

[10567](#)

### UniProt ID

[Q9UI14](#)