

Product Information

MemDX™ Membrane Protein Human GOLT1B (Golgi transport 1B) Full Length

Cat. No.: **MPC3659K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human GOLT1B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GOLT1B

Protein Length

Full length

Protein Class

Transporter

TMD

4

Sequence

MISLTDTQKIGMGLTGFGVFFLFFGMILFFDKALLAIGNVLFVAGLAFVI
GLERTFRFFFQKHMKATGFFLGGVFVVLIGWPLIGMIFEIYGFFLLFRG
FFPVVVGFIIRVPVLGSLNLPGIRSFVDKVGESNNMV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

GOLT1B

Full Name

Golgi transport 1B

Alternative Names

GOLT1B; GCT2; GOT1; GOT1B; CGI-141; YMR292W; vesicle transport protein GOT1B; germ cell tumor 2; hGOT1a; putative NF-kappa-B-activating protein 470; Golgi transport 1B

Gene ID

[51026](#)

UniProt ID

[Q9Y3E0](#)