

Product Information

MemDX™ Membrane Protein Human FAM189A1 (Family with sequence similarity 189 member A1) Full Length

Cat. No.: **MPC3664K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human FAM189A1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

FAM189A1

Protein Length

Full length

Protein Class

Transporter

TMD

4

Sequence

MPPAGGPRAPRPAALPRSLSRLRECPGRSRIVLALGATQMALGCLIVAVS
FAALALTTSARVRHSCPFWAGFSVLLSGLIGVVS WKRPLSLVITFFMLLS
AVCVMLNLAGSILSCQNAQLVNSLEGCQLIKFDSVEVCVCCELQHQS SGC
SNLGETLKLNLPLQENCNAVRLTLKDLLFSVCALNVLSTIVCALATAMCCM
QMVS SDVLQMFLPQRSH PANPTCVTPHGTVLHQTLDFDEFIPPLPPPPYY
PPEYTCTPSTEAQRLHLDFA PSPFGTLYDVAINSPGLLYPAELPPPYEA
VVGQPPASQVTSIGQQVAESSSGDPNTSAGFSTPVPADSTSLLVSEGTAT
PGSSPSDPDGPVGAPAPSEPALPPGHVSPEDPGMG SQVQPGPGRVSRSTSD
PTLCTSSMAGDASSHRPSCSQDLEAGLSEAVPGSASMSRSATAACRAQLS
PAGDPDTWKTDQRPTPEPFPATSKERPRSLVDSKAYADARVLVAKFLEHS
HCALPTEAQHMVGAMRLAVTNEERLEEEAVFGADVLDQV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

FAM189A1

Full Name

Family with sequence similarity 189 member A1

Alternative Names

FAM189A1; TMEM228; protein FAM189A1; transmembrane protein 228; Family with sequence similarity 189 member A1

Gene ID

[23359](#)

UniProt ID

[O60320](#)