

Product Information

MemDX™ Membrane Protein Human GPR137C (G protein-coupled receptor 137C) Full

Length

Cat. No.: **MPC3667K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human GPR137C membrane protein expressed in Baculovirus/Insect expression system. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GPR137C

Protein Length

Full length

Protein Class

GPCR

TMD

7

Sequence

MRVSVPGPAAAAAPAAGREPSTPGGGSGGGGAVAAASGAAVPGSVQLALS
VLHALLYAALFAFAYLQLWRLLLYRERRLSYQSLCLFLCLLWAALRTTLF
SAAFSLSGSLPLLPPAHLHFFPHWLLYCFPSCLQFSTLCLLNLYLAEVI
CKVRCATELDRHKILLHLGFIMASLLFLVNLTCAMLVHGDVPENQLKWT
VFVRALINDSLFILCAISLVCYICKITKMSSANVYLESKGMSLCQTVVVG
SVWILLYSSRACYNLVVVTISQDTLESPFNYGWDNLSDKAHVEDISGEEY
IVFGMVLFLWEHVPAWSVVLFFRAQRLNQNLAPAGMINSHSYSSRAYFFD
NPRRYDSDDLPRLGSSREGSLPNSQSLGWYGMTGCGSSSYTVTPHLNG
PMTDTAPLLFTCSNLDLNNHSLYVTPQN

Product Description

Expression Systems

Baculovirus/Insect expression system

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GPR137C

Full Name

G protein-coupled receptor 137C

Alternative Names

GPR137C; TM7SF1L2; integral membrane protein GPR137C; G protein-coupled receptor TM7SF1L2; transmembrane 7 superfamily member 1-like 2 protein; G protein-coupled receptor 137C

Gene ID

[283554](#)

UniProt ID

[Q8N3F9](#)