

Product Information

MemDX™ Membrane Protein Human STX17 (Syntaxin 17) Full Length

Cat. No.: **MPC3712K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human STX17 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

STX17

Protein Length

Full length

Protein Class

Transporter

TMD

2

Sequence

MSEDEEKVKLRRLEPAIQKFIVPTDLERLRKHQINIEKYQRCRIWDK
LHEEHINAGRTVQQLRNIREIEKLCLKVRKDDLVLKRMIDPVKEEASA
ATAEFLQLHLESVEELKKQFNDEETLLQPPLTRSMTVGGAFHTTEAEASS
QSLTQIYALPEIPQDQNAAESWETLEADLIELSQLVTDLSLLVNSQQEKI
DSIADHVNSAAVNVEEGTKNLGKAAYKLAALPVAGALIGGMVGGPIGLL
AGFKVAGIAAALGGVVGFTGGKLIQRKKQKMMMEKLTSSCPDLPSQTDKK
CS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

STX17

Full Name

Syntaxin 17

Alternative Names

STX17; syntaxin-17; Syntaxin 17

Gene ID

[55014](#)

UniProt ID

[P56962](#)