

# Product Information

## MemDX™ Membrane Protein Human LRRC59 (Leucine rich repeat containing 59) Full

### Length

Cat. No.: **MPC3715K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human LRRC59 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

LRRC59

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MTKAGSKGGNLRDKLDGNELDSLSDLNEVPVKELAALPKATILDLSCHK  
LTTLPSPDFCGLTHLVKLDLSKNKLQQLPADFGRLVNLQHLDLLNNKLVTL  
PVSFAQLKNLKWLDLKDNPVLAKEVAGDCLDEKQCKQCANKVLQHMKA  
VQADQERERQRRLEVEREAEEKKREKQRAKEAQERELRKREKAEKERRR  
KEYDALKAAREQEKKPKKEANQAPKSKSGSRPRKPPPRKHTRSWAVLKL  
LLLLLLFGVAGGLVACRVTELQQQLCTSVNTIYDNAVQGLRRHEILQWV  
LQTDSQQ

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

LRRC59

**Full Name**

Leucine rich repeat containing 59

**Alternative Names**

LRRC59; p34; PRO1855; leucine-rich repeat-containing protein 59; ribosome-binding protein p34; Leucine rich repeat containing 59

**Gene ID**

[55379](#)

**UniProt ID**

[Q96AG4](#)