

Product Information

MemDX™ Membrane Protein Human WBP1L (WW domain binding protein 1 like) Full Length

Cat. No.: **MPC3722K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human WBP1L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

WBP1L

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MPFLLGLRQDKEACVGTNNQSYICDTGHCCGQSQCCNYYYELWWFWLVWT
IIILSCCCVCHHRAKHRLQAQQRQHEINLIAYREAHNYSALPFYFRFL
PNYLLPPYEEVVNRPPTPPPYSAFQLQQQQLLPPQCGPAGGSPPGIDPT
RGSQGAQSSPLSEPSRSSTRPPSIADPDPSDLPVDRAATKAPGMEPSGSV
AGLGELDPGAFLDKDAECREELLKDDSSSEHGAPDSKEKTPGRHRRFTGDS
GIEVCVCNRRGHDDDLKEFNTLIDDALDGPLDFCDSCHVRPPGDEEEGLC
QSSEEQAREPGHPLRPPACLLNTINEQDSPNSQSSSSPS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

WBP1L

Full Name

WW domain binding protein 1 like

Alternative Names

WBP1L; OPA1L; OPAL1; C10orf26; outcome predictor in acute leukemia 1; WW domain binding protein 1 like

Gene ID

[54838](#)

UniProt ID

[Q9NX94](#)