

# Product Information

## MemDX™ Membrane Protein Human SHISAL1 (Shisa like 1) Full Length

Cat. No.: **MPC3723K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SHISAL1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

SHISAL1

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MTSCGQQSLNVLAVLFSLLFSAVLSAHFRVCEPYTDHKGRYHFGFHCPR  
SDNKTILCCHHNNTVFKYCCNETEFQAVMQANLTASSEGYMHNNYTALL  
GVWIYGFFVLMLLVLDLLYYSAMNYDICKVYLARWGIQGRWMKQDPRRWG  
NPARAPRPGQRAPQPQPPGPLPQAPQAVHTLRGDAHSPPLMTFQSSSA

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

#### Form

Liquid

#### Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

## Target

### Target Protein

SHISAL1

### Full Name

Shisa like 1

### Alternative Names

SHISAL1; KIAA1644; protein shisa-like-1; Shisa like 1

### Gene ID

[85352](#)

### UniProt ID

[Q3SXP7](#)