

# Product Information

## MemDX™ Membrane Protein Human SHISAL2B (Shisa like 2B) Full Length

Cat. No.: **MPC3724K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SHISAL2B membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

SHISAL2B

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

1

#### Sequence

MSEASRLCSGYYSLNQSFVEPFQCPRRGEGAALQYCCGFADLKYYCCSEPG  
SYFPYKHSYMWLSIGALIGLGIAALVLLAFVISVCVLCYLFLYTKPQRL  
DTGLKLQHLEASSTQEGKSNGKTKALNSNAASNATNETYYEADDIIQEK  
MDATQIHAIY

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

#### Form

Liquid

#### Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

## Target

### Target Protein

SHISAL2B

### Full Name

Shisa like 2B

### Alternative Names

SHISAL2B; FAM159B; protein shisa-like-2B; family with sequence similarity 159 member B; membrane protein  
FAM159B; Shisa like 2B

### Gene ID

[100132916](#)

### UniProt ID

[A6NKW6](#)