

# Product Information

## MemDX™ Membrane Protein Human SIGLEC9 (Sialic acid binding Ig like lectin 9) Full

### Length

Cat. No.: **MPC3750K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SIGLEC9 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

SIGLEC9

#### Protein Length

Full length

#### Protein Class

Cell adhesion

#### TMD

1

#### Sequence

MLLLLLPLLLWGRERAEGQTSKLLTMQSSVTVQEGLCVHVPCSFSYPSHGW  
IYPGPVVHGYWFRGANTDQDAPVATNNPARAVWEETRDRFHLLGDPHTK  
NCTLSIRDARRSDAGRYFFRMEKGSIKWNYKHHRLSVNVTALTHRPNILI  
PGTLESGCPQNLTCVWPWACEQGTTPMISWIGTSVSPDPSTTRSSVLT  
IPQPQDHGTSLTQVTFPGASVTTNKTVHLNVSYPPQNLMTVFQGDGTV  
STVLGNGSSLSLPEGQSLRLVCAVDAVDSNPPARLSLSWRGLTLCPSQPS  
NPGVLELPWVHLRDAAEFTCRAQNPLGSQQVYLVSLQSKATSGVTQGVV  
GGAGATALVFLSFCVIFVVVRSCRKKSARPAAGVGDTGIEDANAVRGAS  
QGPLETPWAEDSPDPQPPASARSSVGEGELQYASLSFQMVKPWDSRGQE  
ATDTEYSEIKIHR

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

**Protein Format**

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

**Form**

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

SIGLEC9

**Full Name**

Sialic acid binding Ig like lectin 9

**Alternative Names**

SIGLEC9; CD329; CDw329; FOAP-9; siglec-9; OBBP-LIKE; sialic acid-binding Ig-like lectin 9; protein FOAP-9; Sialic acid binding Ig like lectin 9

**Gene ID**

[27180](#)

**UniProt ID**

[Q9Y336](#)