

Product Information

MemDX™ Membrane Protein Human ZDHHC20 (Zinc finger DHHC-type palmitoyltransferase

20) Full Length

Cat. No.: **MPC3767K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ZDHHC20 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ZDHHC20

Protein Length

Full length

Protein Class

Transferase

TMD

4

Sequence

MAPWTLWRCCQRVVGWVPVLFITFVVWSYYAYVVELCVFTIFGNEENGK
TVVYLVAFHLFFVMFVWSYWMTIFTSPASPSKEFYLSNSEKERYEKEFSQ
ERQQEILRRAARALPIYTTSASKTIRYCEKCQLIKPDRAHHCSACDSCIL
KMDHHCPWVNNCVGFSSNYKFFLLFLLYSLLYCLFVAATVLEYFIKFWTNE
LTDTRAKFHVLFVFFVSAMFFISVLSLFSYHCWLVGKNRTTIESFRAPTF
SYGPDGNGFSLGCSKNWRQVFGDEKKYWLLPIFSSSLGDGCSFPTRLVGMD
PEQASVTNQNEYARSSGSNQPFPIKPLSESKNRLLDSESQWLENGAEEGI
VKSGTNNHVTVAIEN

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ZDHHC20

Full Name

Zinc finger DHHC-type palmitoyltransferase 20

Alternative Names

ZDHHC20; DHHC20; DHHC-20; 4933421L13Rik; palmitoyltransferase ZDHHC20; DHHC domain-containing cysteine-rich protein 20; acyltransferase ZDHHC20; probable palmitoyltransferase ZDHHC20; zinc finger DHHC domain-containing protein 20; zinc finger DHHC-type containing 20; Zinc finger DHHC-type palmitoyltransferase 20

Gene ID

[253832](#)

UniProt ID

[Q5W0Z9](#)