

Product Information

MemDX™ Membrane Protein Human DPY19L1 (Dpy-19 like C-mannosyltransferase 1) Full Length

Cat. No.: **MPC3772K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human DPY19L1 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

DPY19L1

Protein Length

Full length

Protein Class

Transferase

TMD

11

Sequence

MEGRPPPEGRPPPRPRTGRAPRGRRRAVFAAVLHWSHITHLFENDRHFSH
LSTLEREMAFRTEMGLYYSYFKTIVEAPSFLNGVWMIMNDKLTEYPLVIN
TLKRFNLYPEVILASWYRIYTKIMDLIGIQTKICWTVTRGEGLSPIESCE
GLGDPACFYVAVIFILNGLMMALFFIYGTYLSGSRLGGLVTVLCFFFNHG
ECTRVMWTPPLRESFSYPFLVLQMLLVTHILRATKLYRGSLIALCISNVF
FMLPWQFAQFVLLTQIASLFAVYVVGIDICKLRKIIYIHMISLALCFVL
MFGNSMLLTSYYASSLVIIWGILAMKPHFLKINVSESLWVIQGCFWLFG
TVILKYLTSKIFGIADDAHIGNLLTSKFFSYKDFDTLLYTCAAEDFMEK
ETPLRYTKTLLLPPVVLVVFVAIVRKIISDMWGVLAQQTTHVRKHQFDHGE
LVYHALQLLAYTALGILIMRLKFLTPHMCVMASLICSRLFGWLFCKVH
PGAIVFAILAAMSIQGSANLQTQWNIVGEFSNLPQEELIEWIKYSTKPD
VFAGAMPTMASVKLSALRPVNHPHYEDAGLRARTKIVYSMYSRKAEEV
KRELKLVNYYILEESWCVRRSKPGCSMPFIWDVEDPANAGKTPLCNLL
VKDSKPHFTTVFQNSVYKVLEVVE

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

DPY19L1

Full Name

Dpy-19 like C-mannosyltransferase 1

Alternative Names

DPY19L1; probable C-mannosyltransferase DPY19L1; DPY-19-like protein 1; protein dpy-19 homolog 1; Dpy-19 like C-mannosyltransferase 1

Gene ID

[23333](#)

UniProt ID

[Q2PZI1](#)