

Product Information

MemDX™ Membrane Protein Human LDLRAD4 (Low density lipoprotein receptor class A domain containing 4) Full Length

Cat. No.: **MPC3774K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human LDLRAD4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

LDLRAD4

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MPEAGFQATNAFTECKFTCTSGKCLYLGSLVCNQQNDCGDNDSDEENCLLV
TEHPPPGIFNSELEFAQIIIVVVVTVMVVIVCLLNHYKVSTRSFINRP
NQSRRREDGLPQEGCLWPSDSAAPRLGASEIMHAPRSRDRFTAPSFQIRD
RFSRFQPTYPYVQHEIDLPTISLSDGEEPPPYQGPCTLQLRDPEQQMEL
NRESVRAPPNRTIFDSDLIDIAMYSGGPCPPSSNSGISASTCSSNGRMEG
PPPTYSEVMGHHHPGASFLHHQRSNAHRGSRLQFQQNNAESTIVPIKGKDR
KPGNLV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

LDLRAD4

Full Name

Low density lipoprotein receptor class A domain containing 4

Alternative Names

LDLRAD4; C18orf1; low-density lipoprotein receptor class A domain-containing protein 4; clone 22; Low density lipoprotein receptor class A domain containing 4

Gene ID

[753](#)

UniProt ID

[O15165](#)