

Product Information

MemDX™ Membrane Protein Human SIGLEC14 (Sialic acid binding Ig like lectin 14) Full Length

Cat. No.: **MPC3777K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SIGLEC14 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SIGLEC14

Protein Length

Full length

Protein Class

Cell adhesion

TMD

1

Sequence

MLPLLLLPLLWGGSLQEKPVYELQVQKSVTVQEGLCVLVPCSFSPWRSW
YSSPPLYVYWFRDGEIPYYAEVVATNNPDRRVKQGRFRLLGDVQKKN
CSLSIGDARMEDTGSYFFRVERGRDVKYSYQQNKLNLVETALIEKPDHF
LEPLESGRPTLRLCSLPGSCEAGPPLTFSTGNALSPLDPETTRRSELT
TPRPEDHGTNLTCQVKRQGAQVTTERTVQLNVSYAPQNLAISSFRNGTG
TALRILSNGMSVPIQEGQSLFLACTVDSNPPASLSWFREGKALNPSQTSM
SGTLELPNIGAREGGFTCRVQHPLGSQHLSFILSVQRSSSSCICVTEKQ
QGSWPLVLTIRGALMGAGFLLTYGLTWIYYTRCGGPQQSRAERPG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SIGLEC14

Full Name

Sialic acid binding Ig like lectin 14

Alternative Names

SIGLEC14; sialic acid-binding Ig-like lectin 14; siglec-14; Sialic acid binding Ig like lectin 14

Gene ID

[100049587](#)

UniProt ID

[Q08ET2](#)