

Product Information

MemDX™ Membrane Protein Human GLT8D2 (Glycosyltransferase 8 domain containing 2)

Full Length

Cat. No.: **MPC3813K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human GLT8D2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GLT8D2

Protein Length

Full length

Protein Class

Transferase

TMD

1

Sequence

MALLRKINQVLLFLLIVTLCVILYKKVHKGTVPKNDADDESETPEEEEEE
IPVVICAAAGRMGATMAAINSISNTDANILFYVVGLRNTLTRIRKWIEH
SKLREINFKIVEFNPMVLKGKIRPDSSRPELLQPLNFRFYLPPLIHQHE
KVIYLDLDDVIVQGDIQELYDTTLALGHAAAFSDDCDLPQAQDINRLVGLQ
NTYMGYLDYRKKAIKDLGISPSTCSFNPGVIVANMTEWKHQIRITKQLEKW
MQKNVEENLYSSSLGGGVATSPMLIVFHGKYSTINPLWHIRHLGWNPDPAR
YSEHFLQEAKLLHWNGRHKPWDFPSVHNDLWESWFPDPAGIFKLNHHS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GLT8D2

Full Name

Glycosyltransferase 8 domain containing 2

Alternative Names

GLT8D2; glycosyltransferase 8 domain-containing protein 2; Glycosyltransferase 8 domain containing 2

Gene ID

[83468](#)

UniProt ID

[Q9H1C3](#)