

Product Information

MemDX™ Membrane Protein Human MPZL3 (Myelin protein zero like 3) Full Length

Cat. No.: **MPC3817K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human MPZL3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

MPZL3

Protein Length

Full length

Protein Class

Cell adhesion

TMD

1

Sequence

MQQRGAAGSRGCALFPLLGVLFQGVYIVFSLEIRADAHVRGYVGEKIKL
KCTFKSTSDVTDKLTIDWTYRPPSSSHTVSIFHYQS FQYPTTAGTFRDRI
SWVGNVYKGDASISNPTIKDNGTFSCAVKNPPDVHHNIPMTELTVTER
GFGTMLSSVALLSILVFVPSAVVVALLLVRMGRKAAGLKKRSRSGYKKSS
IEVSDDTDQEEEEACMARLCVRCAECLDSYEETY

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

MPZL3

Full Name

Myelin protein zero like 3

Alternative Names

MPZL3; myelin protein zero-like protein 3; Myelin protein zero like 3

Gene ID

[196264](#)

UniProt ID

[Q6UWV2](#)