

Product Information

MemDX™ Membrane Protein Human CEACAM19 (CEA cell adhesion molecule 19) Full

Length

Cat. No.: **MPC3846K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CEACAM19 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CEACAM19

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEIPMG TQGCF SKSLLSAS ILVLWMLQGSQAALYIQKIPEQPQKNQDLL
LSVQGV PDTFQDFN WYLGEETYGGTRLFTYIPGIQRPQRDGSAMGQRDIV
GFPNGS MLLRRAQPTDS GTYQVAITINSEWTMKAKTEVQVAEKNKELPST
HLPTNAG ILAATIIGSLAAGALLISCIAYLLVTRNWRGQSHRLPAPRGQG
SLSILCSAVSPVPSVTPSTWMATTEKPELGPAHDAGDNNIYEVMPSPVLL
VSPISDTRSINPARPLPTPPHLQAEPENHQYQQDLLNPD PAPYCQLVPTS

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

CEACAM19

Full Name

CEA cell adhesion molecule 19

Alternative Names

CEACAM19; CEAL1; CEACM19; carcinoembryonic antigen-related cell adhesion molecule 19; carcinoembryonic antigen related cell adhesion molecule 19; carcinoembryonic antigen-like 1; CEA cell adhesion molecule 19

Gene ID

[56971](#)

UniProt ID

[Q7Z692](#)