

Product Information

MemDX™ Membrane Protein Human CD99L2 (CD99 molecule like 2) Full Length

Cat. No.: **MPC3865K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CD99L2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CD99L2

Protein Length

Full length

Protein Class

Cell adhesion

TMD

1

Sequence

MVAWRS AFLVCLAFSLATLVQRGSGDFDDFNLEDAVKETSSVKQPWDHTT
TTTTNRPGTTTRAPAKPPGSGLDLADALDDQDDGRRKPGIGGRERWNHVTT
TTKRPVTTTRAPANTLGNDFDLADALDDRNRDDGRRKPIAGGGGFSDKDL
EDIVGGGEYKPKDGKGDGRYGSNDDPGSGMVAEPGTIAGVASALAMALIG
AVSSYISYQQKKFCFSIQQLNADYVKGENLEAVVCEEPQVKYSTLHTQS
AEPPPPPEPARI

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

CD99L2

Full Name

CD99 molecule like 2

Introduction

This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants.

Alternative Names

CD99L2; CD99B; MIC2L1; CD99 antigen-like protein 2; MIC2 like 1; MIC2-like protein 1; CD99 molecule like 2

Gene ID

[83692](#)

UniProt ID

[Q8TCZ2](#)