

Product Information

MemDX™ Membrane Protein Human SUS4 (Sushi domain containing 4) Full Length

Cat. No.: **MPC3908K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SUS4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SUS4

Protein Length

Full length

Protein Class

Immunity

TMD

1

Sequence

MYHGMNPSNGDGFLEQQQQQQQPQSPQRLLAVILWFQLALCFGPAQLTGG
FDDLQVCADPGIPENGFRTPSGGVFFEGSVARFHCQDGFKLKGATKRLCL
KHFNGTLGWIPSDNSICVQEDCRIPQIEDAEIHNKTYRHGEKLIITCHEG
FKIRYPDLHNMVSLCRDDGTWNNLPICQGCLRPASSNGYVNISELQTSF
PVGTVISYRCFPFGKLDGSAYLECLQNLIWSSPPRCLALEVCPLPPMVS
HGDFVCHPRPCERYNHGTVEFYCDPGYSLTSDYKYITCQYGEWFPSYQV
YCIKSEQTWPSTHETLLTTWKIVAFTATSVLLVLLLVLARMFQTKFAH
FPPRGPPRSSSDPDFVVVDGVPVMLPSYDEAVSGGLSALGPGYMASVGQ
GCPLPVDDQSPPAYPGSGD TDTGPGESETCDSVSGSSELLQSLYSPPRCQ
ESTHPASDNPDIIASTAEVASTSPGIDIADEIPLMEEDP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SUSD4

Full Name

Sushi domain containing 4

Alternative Names

SUSD4; PRO222; sushi domain-containing protein 4; YHGM196; testicular secretory protein Li 55; Sushi domain containing 4

Gene ID

[55061](#)

UniProt ID

[Q5VX71](#)