

Product Information

MemDX™ Membrane Protein Human GSG1L (GSG1 like) Full Length

Cat. No.: **MPC3917K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human GSG1L membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

GSG1L

Protein Length

Full length

Protein Class

Receptor

TMD

4

Sequence

MKTSRRGRALLAVALNLLALLFATTAFLTTHWCQGTQRVPKPGCGQGGRA
NCPNSGANATANGTAAPAAAAAATASNGPPGGALYSWETGDDRFLFRN
FHTGIWYSCEEELSGLGEKCRSFIDLAPASEKGVWLSVVSEVLYILLV
VGFSLMCLELFHSSNVIDGLKLNFAAAVFTVLSGLLGMAHMMYTQVFQV
TVSLGPEDWRPHSWDYGWSFCLAWGSFTCCMAASVTTLNSYTKTVIEFRH
KRKVFEQGYREEPTFIDPEAIKYFRERMEKRDGSEEDFHLD CRHERYPAR
HQPHMADSWPRSSAQEAPELNRQCWVLGHVV

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

GSG1L

Full Name

GSG1 like

Alternative Names

GSG1L; PRO19651; germ cell-specific gene 1-like protein; GSG1-like protein; KTSR5831; GSG1 like

Gene ID

[146395](#)

UniProt ID

[Q6UXU4](#)