

Product Information

MemDX™ Membrane Protein Human CDHR4 (Cadherin related family member 4) Full Length

Cat. No.: **MPC3934K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human CDHR4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

CDHR4

Protein Length

Full length

Protein Class

Cell adhesion

TMD

1

Sequence

MVLLRLLVFLFAPVVSDLCSLPCFINVSESQGPGLVQLFLSFNCSSYTPT
PTLELLNVQPPTTFFNPPSLARWQGTYYGKLTLSAQLDALMVNHYKVQ
LKFTCGNHVMEGSLSDVDVQRDLSHIQCAGQFASPAGEMIQVPETVTPGAR
LYTLLLPGLELHGAQMSIISAQDLPHFPGPFSINEQGWLQAPSQGLLGQA
QKVFQLQISVSFGQRQSCQGMVIVKVLPVPSSQVSFLEQAQNTIPENLA
PGSEVVQVQARGVDLRYEILSPVPSPLFSIGRADGVVRTTTPELARTSG
TAVSRLQVKAFEQQLWASAKLNLTMNVQLVNLWPPRCLPALLVSQIPET
APVGTVLNLTLCEDPDSVGATLDYKLWFRSSNPASLCYDRVLEVNLATL
DCDTPGACFQHAASILVLDGGQPQMTTEVPVLMVTPINEFSPACAPRTF
RVQEDAAPHTLLGSVVGTDMDYPHDNIEYYTSGGPTTFAVDRLSGEVHLL
GPLDYEQQRLYRLTVLVIDHGQDQNPNNHLSGSCITIEVEDVNDHAPEC
EPPFQELTIYAPLGRSVEVTKMSCIQEPQRLIYSYSIVGGNSQNRFIL
QGAILVHSDLVLGPFWPEQPRTYELLICVADAGPSTPHLSTTATIIVHLV
PRRASTVATSTHRTTVPSTMTMLVTDTEAFWQPQPWFVVLTATGALLL
LALGWLLGRLLQGLAQLLQAPSKPAQALLNSIQGTEGSIEGFLEAPKME
MSQAPSSVMSLHFDGRAQDSRTGRDYLFNTHTGARRWL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

CDHR4

Full Name

Cadherin related family member 4

Alternative Names

CDHR4; CDH29; PRO34300; cadherin-related family member 4; Cadherin-like protein UNQ9392/PRO34300; VLLR9392; cadherin-like 29; cadherin-like protein 29; Cadherin related family member 4

Gene ID

[389118](#)

UniProt ID

[A6H8M9](#)