

Product Information

MemDX™ Membrane Protein Human TOR4A (Torsin family 4 member A) Full Length

Cat. No.: **MPC3983K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human TOR4A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

TOR4A

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MDRGQPSLEPAAAAPRASGRCVIAPVRAVLRLRRRVCVLRKRRLQLPGGG
PDVGTGAPRPGCSPRAPRADLDQPKFFTFDSPAELPSRTPRKRRRSRLV
LYPETSRYRPRVEHRSRAQRCLLLLVAIVGFQVLNAIENLDDNAQRYDL
DGLEKALQRAVFGQPAAVSRIVALMRDYLATHVHSRPLLLALHGPGSGVGK
SHVGRLLARHFRSVLEDSALVLQYHARHHCPAARAAQDCREELARRVADV
VARAEAEKTPLLVLDDVELMPRPLLDLHGFLQPQRSHHFHNAIYVLLS
GAGGAEVTRFVLQNASRALPLRPDGFSAEAAAAQAEEDLRASLLAVLSR
EHPLWQAAAIVPFLLLDKRDVVSCFRDEMAGEGFFPDQARAENLAAQLSF
YRVAGREFAVTGCKQVVATVNLL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

TOR4A

Full Name

Torsin family 4 member A

Alternative Names

TOR4A; C9orf167; torsin-4A; torsin family protein C9orf167; Torsin family 4 member A

Gene ID

[54863](#)

UniProt ID

[Q9NXH8](#)