

Product Information

MemDX™ Membrane Protein Human STX4 (Syntaxin 4) Full Length

Cat. No.: **MPC4023K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human STX4 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

STX4

Protein Length

Full length

Protein Class

Transporter

TMD

1

Sequence

MRDRTHELRQGDDSSDEEDKERVALLVHPGTARLGSPDEEFFHKVRTIRQ
TIVKLGKQVQLEKQQVTILATPLPEESMKQELQNLRLDEIKQLGREIRLQ
LKAIEPQKEEADENYSVNTRMRKTQHGVLSQQFVELINKCNSMQSEYRE
KNVERIRRLKITNAGMVSDEELEQMLDSGQSEVFVSNILKDTQVTRQAL
NEISARHSEIQQLERSIRELHDIFTFLATEVEMQGEMINRIEKNILSSAD
YVERGQEHVKTALNQKKARKKKVLIAICVSITVLLAVIIGVTVVG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

STX4

Full Name

Syntaxin 4

Alternative Names

STX4; STX4A; p35-2; syntaxin-4; renal carcinoma antigen NY-REN-31; syntaxin 4A (placental); Syntaxin 4

Gene ID

[6810](#)

UniProt ID

[Q12846](#)