

Product Information

MemDX™ Membrane Protein Human SIGLEC15 (Sialic acid binding Ig like lectin 15) Full Length

Cat. No.: **MPC4167K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human SIGLEC15 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

SIGLEC15

Protein Length

Full length

Protein Class

Receptor

TMD

1

Sequence

MEKSIWLLACLAWLPTGFSFVRTKIDTTENLLNTEVHSSPAQRWSMQVPP
EVSAEAGDAAVLPCTFTHPHRHYDGPLTAIWRAGEPYAGPQVFRCAAARG
SELCQTALSLHGRFRLGNPRRNDLSLRVERLALADDRRYFCRVEFAGDV
HDRYESRHHGVRLHVTAAAPRIVNISVLPSPAHAFFRALCTAEGEPPPALAWS
GPALGNSLAAVRSPREGHGHVTAELPALTHDGRYTCTAANSLGRSEASV
YLFRFHGASGASTVALLLGALGFKALLLLGVLAARAARRRPEHLDTPDTP
PRSQAQESNYENLSQMNPRSPPATMCSP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

SIGLEC15

Full Name

Sialic acid binding Ig like lectin 15

Alternative Names

SIGLEC15; CD33L3; HsT1361; SIGLEC-15; CD33 antigen-like 3; CD33 molecule-like 3; Sialic acid binding Ig like lectin 15

Gene ID

[284266](#)

UniProt ID

[Q6ZMC9](#)