

# Product Information

## MemDX™ Membrane Protein Human YIPF6 (Yip1 domain family member 6) Full Length

Cat. No.: **MPC4168K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human YIPF6 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

### Product Specifications

#### Host Species

Human

#### Target Protein

YIPF6

#### Protein Length

Full length

#### Protein Class

Receptor

#### TMD

5

#### Sequence

MAEAEESPGDPTASPRPLFAGLSDISISQDIPVEGEITIPMRSRIREFD  
SSTLNESVRNTIMRDLKAVGKKFMHVLYPRKSNTLLRDWDLWGPLILCVT  
LALMLQRDSADSEKDGGPQFAEVFVIVWFGAVTITLNSKLLGGNISFFQS  
LCVLGYCILPLTVAMLICRLVLLADPGPVNFMVRLFFVIVMFAWSIVAST  
AFLADSQPPNRRALAVYPVFLFYFVISWMILTFTPQ

### Product Description

#### Expression Systems

HEK293

#### Tag

Based on specific requirements

#### Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

#### Form

Liquid

**Storage**

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

**Target****Target Protein**

YIPF6

**Full Name**

Yip1 domain family member 6

**Alternative Names**

YIPF6; FinGER6; protein YIPF6; YIP1 family member 6; Yip1 domain family member 6

**Gene ID**

[286451](#)

**UniProt ID**

[Q96EC8](#)