

Product Information

MemDX™ Membrane Protein Human EXTL2 (Exostosin like glycosyltransferase 2) Full

Length

Cat. No.: **MPC4195K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human EXTL2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

EXTL2

Protein Length

Full length

Protein Class

Transferase

TMD

1

Sequence

MRCCHICKLPGRVMGIRVLRSLVVLVLLL VAGALTALLPSVKEDKMLM
LRREIKSQGKSTMDSFTLMQTYNRTDLLL KLLNHYQAVPNLHKVIVVWN
NIGEKAPDELWNSLGPHPIPVIFKQQTANRMRNRLQVFPELETNAVLMVD
DDTLISTPDLVFAFSVWQQFPDQIVGFVPRKHVSTSSGIYSYGSFEMQAP
GSGNGDQYSMVLIGASFFNSKYLELFQRQPAAVHALIDDTQNCDDIAMNF
IIAKHIGKTS GIFVKPVNMDNLEKETNSGYSGMWHRAEHALQRSYCINKL
VNIYDSMPLRYSNIMISQFGFPYANYKRKI

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

EXTL2

Full Name

Exostosin like glycosyltransferase 2

Alternative Names

EXTL2; EXTR2; exostosin-like 2; EXT-related protein 2; alpha-1,4-N-acetylhexosaminyltransferase EXTL2; alpha-GalNAcT EXTL2; exostosin-like 2; glucuronyl-galactosyl-proteoglycan 4-alpha-N-acetylglucosaminyltransferase; processed exostosin-like 2; Exostosin like glycosyltransferase 2

Gene ID

[2135](#)

UniProt ID

[Q9UBQ6](#)