

Product Information

MemDX™ Membrane Protein Human HIGD1A (HIG1 hypoxia inducible domain family member 1A) Full Length

Cat. No.: **MPC4313K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human HIGD1A membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

HIGD1A

Protein Length

Full length

Protein Class

Transporter

TMD

2

Sequence

MSTDTGVSLPSYEEDQGSKLIRKAKEAPFVPVGIAGFAAIVAYGLYKLKS
RGNTKMSIHLIHRVAAQGFVVGAMTVGMGYSMYREFWAKPKP

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target

Target Protein

HIGD1A

Full Name

HIG1 hypoxia inducible domain family member 1A

Alternative Names

HIGD1A; HIG1; RCF1a; HIG1 domain family member 1A, mitochondrial; HIG1 domain family, member 1A; RCF1 homolog A; hypoxia inducible gene 1; hypoxia-inducible gene 1 protein; HIG1 hypoxia inducible domain family member 1A

Gene ID

[25994](#)

UniProt ID

[Q9Y241](#)