

Product Information

MemDX™ Membrane Protein Human NDFIP2 (Nedd4 family interacting protein 2) Full

Length

Cat. No.: **MPC4325K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human NDFIP2 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

NDFIP2

Protein Length

Full length

Protein Class

Receptor

TMD

3

Sequence

MARRRSQRVCASGPSMLNSARGAPELLRGTATNAEVSAAAAGATGSEELP
PGDRGCRNGGGRGPAATTSSTGVAVGAEHGEDSLSRKPDPEPGRMDHHQP
GTGRYQVLLNEEDNSESSAIEQPPTSNPAPQIVQAASSAPALETDSPPP
YSSITVEVPTTSDTEVYGEFYVPPPYSVATSLPTYDEAEKAKAAAMAAA
AAETSQRIQEEECPPRDDFSDADQLRVGNDGIFMLAFFMAFIFNWLGFCL
SFCITNTIAGRYGAICGFGLSLIKWILIVRFSDYFTGYFNGQYWLWWIFL
VLGLLLFFRGFVNYLKVRNMSESMAAAHRTRYFFLL

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

NDFIP2

Full Name

Nedd4 family interacting protein 2

Alternative Names

NDFIP2; N4WBP5A; MAPK-activating protein PM04 PM05 PM06 PM07; NEDD4 WW domain-binding protein 5A; NF-kappa-B-activating protein 413; putative MAPK-activating protein PM04/PM05/PM06/PM07; putative NF-kappa-B-activating protein 413; Nedd4 family interacting protein 2

Gene ID

[54602](#)

UniProt ID

[Q9NV92](#)