

Product Information

MemDX™ Membrane Protein Human ERGIC3 (ERGIC and golgi 3) Full Length

Cat. No.: **MPC4340K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human ERGIC3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

ERGIC3

Protein Length

Full length

Protein Class

Transporter

TMD

2

Sequence

MEALGKLKQFDAYPKTLEDFRVKTCGGATVTIVSGLLMLLLFLSELQYYL
TTEVHPELYYDKSRGDKLKINIDVLFPHMPCAYLSIDAMDVAGEQQLDVE
HNLFKQRLDKDGIPVSSEAERHELKGVEVTVPDPSLDPDRCESCYGAEA
EDIKCCNTCEDVREAYRRRGWAFKNPDTIEQCRREGFSQKMQEKNEGCQ
VYGFLEVNVKAGNFHFAPGKSFQQSHVHVHDLQSFGLDNINMTHYIQHLS
FGEDYPGIVNPLDHTNVTAPQASMMFQYFVKVPTVYMKVDGEVLRTNQF
SVTRHEKVANGLLDQGLPGVFLYELSPMMVKLTEKHRSFTHFLTGVCA
IIGGMFTVAGLIDSLIYHSARAIQKKIDLGKTT

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

ERGIC3

Full Name

ERGIC and golgi 3

Alternative Names

ERGIC3; Erv46; CGI-54; C2orf47; PRO0989; C20orf47; NY-BR-84; SDBCAG84; dJ477O4.2; endoplasmic reticulum-Golgi intermediate compartment protein 3; endoplasmic reticulum-localized protein ERp43; serologically defined breast cancer antigen 84; serologically defined breast cancer antigen NY-BR-84; ERGIC and golgi 3

Gene ID

[51614](#)

UniProt ID

[Q9Y282](#)