

Product Information

MemDX™ Membrane Protein Human YIPF3 (Yip1 domain family member 3) Full Length

Cat. No.: **MPC4375K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human YIPF3 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

YIPF3

Protein Length

Full length

Protein Class

Receptor

TMD

5

Sequence

MATTAAPAGGARNGAGPEWGGFEENIQGGGSAVIDMENMDDTSGSSFEDM
GELHQRLREEEVDADAADAAAAEEEDGEFLGMKGFKGQLSRQVADQMWWQA
GKRQASRAFSLYANIDILRPYFDVEPAQVRSRLLESMIPIKMVNFPQKIA
GELYGPLMLVFTLVAILLHGMKTSDTIIREGTLMGTAIGTCFGYWLGVSS
FIYFLAYLCNAQITMLQMLALLGYGLFGHCIVLFITYNIHLHALFYLFWL
LVGGLSTLRMVAVLVSRVTGPTQRLLLCGTLAALHMLFLLYLHFAYHKVV
EGILDTLEGPNIPPIQVRVPRDIPAMLPAARLPTTVLNATAKAVAVTLQSH

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

YIPF3

Full Name

Yip1 domain family member 3

Alternative Names

YIPF3; KLIP1; FinGER3; C6orf109; dJ337H4.3; protein YIPF3; killer lineage protein 1; natural killer cell-specific antigen KLIP1; Yip1 domain family member 3

Gene ID

[25844](#)

UniProt ID

[Q9GZM5](#)