

Product Information

MemDX™ Membrane Protein Human XKR8 (XK related 8) Full Length

Cat. No.: **MPC4386K**

This product is for research use only and is not intended for diagnostic use.

This product is a made-to-order Human XKR8 membrane protein expressed in HEK293. The protein is for research use only and is not approved for use in humans or in clinical diagnosis.

Product Specifications

Host Species

Human

Target Protein

XKR8

Protein Length

Full length

Protein Class

Receptor

TMD

8

Sequence

MPWSSRGALLRDLVLGVLTAAFLDLGTDLWAAVQYALGGRYLWAAVLV
ALLGLASVALQLFSWLWLRADPAGLHGSQPPRRCLALLHLLQLGYLYRCV
QELRQGILLVWQQEEPSEFDLAYADFLALDISMLRLFETFLETAPQLTLVL
AIMLQSGRAEYYQWVGICTSFLGISWALLDYHRALRTCLPSKPLLGLGSS
VIYFLWNLLLLWPRVLAVALFSALFPSYVALHFLGLWLVLVLLVWLQGT
FMPDPSSEWLYRVTATILYFSWFNVAEGRTRGRAIIHFAFLSDSILLV
ATWVTHSSWLPSGIPLQLWLPVGCFCFFLGLALRLVYYHWLHPSCCWKPD
PDQVDGARSLLSPEGYQLPQNRRMTHLAQKFFPKAKDEAASPVKG

Product Description

Expression Systems

HEK293

Tag

Based on specific requirements

Protein Format

Detergent or based on specific requirements (Detergent, Liposome, Nanodisc, Polymer, VLP)

Form

Liquid

Storage

Aliquot and store at -20°C or lower. For long term storage, we recommend to store at -72°C or lower. Avoid freeze/thaw cycles.

Target**Target Protein**

XKR8

Full Name

XK related 8

Alternative Names

XKR8; XRG8; hXkr8; XK-related protein 8; X Kell blood group precursor-related family, member 8; X-linked Kx blood group related 8; XK, Kell blood group complex subunit-related family, member 8; XK related 8

Gene ID

[55113](#)

UniProt ID

[Q9H6D3](#)